



Product Datasheet

NUDT6 antibody (orb129707)

Description	Rabbit polyclonal to NUDT6.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human NUDT6(NP_009014.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	P53370
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...



Immunofluorescence analysis of HeLa cells...