

Product Datasheet

N-(p-Coumaroyl) serotonin (orb1297274)

Description	N-(p-Coumaroyl) serotonin
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{19}H_{18}N_2O_3$
Purity	98.78%
MW	322.36
CAS Number	68573-24-0
SMILES	<chem>C(CNC(C=CC1=CC=C(O)C=C1)=O)C=2C=3C(NC2)=CC=C(O)C3</chem>
Expiration Date	12 months from date of receipt.