



Product Datasheet

RAD23B antibody (orb129730)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal to RAD23B. |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic peptide of human Rad23B(NP_002865.1). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Concentration | 1 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 60kDa |
| Uniprot ID | P54727 |
| Dilution Range | WB: 1:500-2000, IHC-P: 1:100-200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of MCF-7



DU145



293T