

Product Datasheet

Hypaconine (orb1297336)

Description	Hypaconine
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{39}NO_8$
Purity	98%
MW	469.57
CAS Number	63238-68-6
SMILES	<chem>O(C)[C@@H]1[C@@]23[C@]4([C@](COC)(CN(C)[C@@]2([C@]([C@@H]4OC)([C@]5(O)[C@@]6([C@]3(C[C@@](O)([C@@H]6O)[C@@H](OC)[C@@H]5O)[H])[H])[H])CC1)[H]</chem>
Expiration Date	12 months from date of receipt.