

Product Datasheet

Caulophyllogenin (orb1297430)

Description	Caulophyllogenin
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{30}H_{48}O_5$
Purity	98.56%
MW	488.7
CAS Number	52936-64-8
SMILES	<chem>C(O)(=O)[C@]12[C@](C=3[C@@](C)(C[C@H]1O)[C@@]4(C)[C@](CC3)([C@]5(C)[C@@](CC4)([C@@](CO)(C)[C@H](O)CC5)[H])[H])(CC(C)(C)CC2)[H]</chem>
Expiration Date	12 months from date of receipt.