

Product Datasheet

LASV inhibitor 3.3 (orb1297774)

Description	LASV inhibitor 3.3
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{30}H_{37}N_3O_2$
Purity	99.55%
MW	471.63
CAS Number	554438-52-7
SMILES	<chem>C(NCC(=O)N1CCN(C(C2=CC=CC=C2)C3=CC=CC=C3)CC1)(=O)C45CC6CC(C4)CC(C5)C6</chem>
Expiration Date	12 months from date of receipt.