

Product Datasheet

RORyt inverse agonist 13 (orb1297973)

Description	RORyt inverse agonist 13
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{23}H_{17}Cl_2F_3N_2O_4$
Purity	98%
MW	513.29
CAS Number	2170477-75-3
SMILES	<chem>OC(Cc1ccc(CNC(Nc2cc(Cl)c(-c(cccc3)c3OC(F)(F)F)c(Cl)c2)=O)cc1)=O</chem>
Expiration Date	12 months from date of receipt.