

Product Datasheet

TLR4-IN-C34 (orb1298877)

Description	TLR4-IN-C34
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{17}H_{27}NO_9$
Purity	99.77%
MW	389.4
CAS Number	40592-88-9
SMILES	<chem>O(C(C)=O)[C@H]1[C@H](OC(C)=O)[C@@H](COC(C)=O)O[C@H](OC(C)C)[C@@H]1NC(C)=O</chem>
Expiration Date	12 months from date of receipt.