

Product Datasheet

Tenofovir alafenamide hemifumarate (orb1299952)

Description	Tenofovir alafenamide hemifumarate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{29}N_6O_5 \cdot 1/2C_4H_4O_4$
Purity	98.37%
MW	534.51
CAS Number	1392275-56-7
SMILES	<chem>OC(/C=C/C(O)=O)=O.CC(C)OC([C@H](C)NP(CO[C@H](C)CN1C2=NC=N C(N)=C2N=C1)(OC3=CC=CC=C3)=O)=O.CC(C)OC([C@H](C)NP(CO[C@ H](C)CN4C5=NC=NC(N)=C5N=C4)(OC6=CC=CC=C6)=O)=O</chem>
Expiration Date	12 months from date of receipt.