

Product Datasheet

Fenoldopam mesylate (orb1300422)

Description	Fenoldopam mesylate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{17}H_{20}ClNO_6S$
Purity	99.19%
MW	401.86
CAS Number	67227-57-0
SMILES	<chem>CS(O)(=O)=O.Oc1=CC=C(C=C1)C1CNCCC2=C(Cl)C(O)=C(O)C=C12</chem>
Expiration Date	12 months from date of receipt.