

## Product Datasheet

### BI-78D3 (orb1300469)

<b>Description</b>	BI-78D3
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{13}H_9N_5O_5S_2$
<b>Purity</b>	99.95%
<b>MW</b>	379.37
<b>CAS Number</b>	883065-90-5
<b>SMILES</b>	<chem>S(C=1N(C=2C=C3C(=CC2)OCCO3)C(=O)NN1)C=4SC(N(=O)=O)=CN4</chem>
<b>Expiration Date</b>	12 months from date of receipt.