

Product Datasheet

TIC10 Isomer (orb1300538)

Description	TIC10 Isomer
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{26}N_4O$
Purity	99.91%
MW	386.49
CAS Number	41276-02-2
SMILES	<chem>C(N1C2=C(C(=O)N3C1=NCC3)CN(CC4=CC=CC=C4)CC2)C5=C(C)C=CC=C5</chem>
Expiration Date	12 months from date of receipt.