

Product Datasheet

Biapenem (orb1300791)

Description	Biapenem
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{15}H_{18}N_4O_4S$
Purity	99.77%
MW	350.39
CAS Number	120410-24-4
SMILES	<chem>C([O-])(=O)C=1N2[C@@]([C@@]([C@@H](C)O)(C2=O)[H])([C@@H](C)C1SC3CN4[N+](C3)=CN=C4)[H]</chem>
Expiration Date	12 months from date of receipt.