

Product Datasheet

Bepotastine Besilate (orb1300797)

Description	Bepotastine Besilate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{27}H_{31}ClN_2O_6S$
Purity	99.86%
MW	547.06
CAS Number	190786-44-8
SMILES	<chem>OS(=O)(=O)c1ccccc1.OC(=O)CCCN1CCC(CC1)O[C@@H](c1ccc(Cl)cc1)c1cccn1</chem>
Expiration Date	12 months from date of receipt.