

Product Datasheet

MGCD-265 analog (orb1300844)

Description	MGCD-265 analog
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{26}H_{20}FN_5O_2S_2$
Purity	99.21%
MW	517.6
CAS Number	875337-44-3
SMILES	<chem>O(C1=C2C(C=C(S2)C3=CN(C)C=N3)=NC=C1)C4=C(F)C=C(NC(NC(CC5=CC=CC=C5)=O)=S)C=C4</chem>
Expiration Date	12 months from date of receipt.