

Product Datasheet

TW-37 (orb1300906)

Description	TW-37
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{33}H_{35}NO_6S$
Purity	98%
MW	573.7
CAS Number	877877-35-5
SMILES	<chem>S(=O)(=O)(C1=C(C(C)(C)C)C=CC=C1)C2=CC=C(NC(=O)C3=CC(CC4=C(C(C)C)C=CC=C4)=C(O)C(O)=C3O)C=C2</chem>
Expiration Date	12 months from date of receipt.