

Product Datasheet

ISRIB (trans-isomer) (orb1301002)

Description	ISRIB (trans-isomer)
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{22}H_{24}Cl_2N_2O_4$
Purity	98.81%
MW	451.34
CAS Number	1597403-47-8
SMILES	<chem>N(C(COC1=CC=C(Cl)C=C1)=O)[C@H]2CC[C@H](NC(COC3=CC=C(Cl)C=C3)=O)CC2</chem>
Expiration Date	12 months from date of receipt.