

## Product Datasheet

### Telotristat Etiprate (orb1301123)

<b>Description</b>	Telotristat Etiprate
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{36}H_{35}ClF_3N_7O_6$
<b>Purity</b>	99.21%
<b>MW</b>	754.15
<b>CAS Number</b>	1137608-69-5
<b>SMILES</b>	<chem>OC(=O)CNC(=O)c1ccccc1.CCOC(=O)[C@@H](N)Cc1ccc(cc1)-c1cc(O[C@H](c2ccc(Cl)cc2-n2ccc(C)n2)C(F)(F)F)nc(N)n1</chem>
<b>Expiration Date</b>	12 months from date of receipt.