

## Product Datasheet

### TG4-155 (orb1301754)

<b>Description</b>	TG4-155
<b>Storage</b>	Store at -20°C for 12 months.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{26}N_2O_4$
<b>Purity</b>	99.14%
<b>MW</b>	394.46
<b>CAS Number</b>	1164462-05-8
<b>SMILES</b>	<chem>C(CNC(/C=C/C1=CC(OC)=C(OC)C(OC)=C1)=O)N2C=3C(C=C2C)=CC=C3</chem>
<b>Expiration Date</b>	12 months from date of receipt.