

Product Datasheet

Treprostinil Sodium (orb1302030)

Description	Treprostinil Sodium
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{23}H_{33}NaO_5$
Purity	99.95%
MW	412.5
CAS Number	289480-64-4
SMILES	<chem>[Na+].CCCCC[C@H](O)CC[C@H]1[C@H](O)C[C@@H]2Cc3c(C[C@H]12)ccc3OCC([O-])=O</chem>
Expiration Date	12 months from date of receipt.