

Product Datasheet

PQ401 hydrochloride (196868-63-0(free base)) (orb1302951)

Description	PQ401 hydrochloride (196868-63-0(free base))
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{18}H_{17}Cl_2N_3O_2$
Purity	99.90%
MW	378.25
SMILES	<chem>Cl.COC1ccc(Cl)cc1NC(=O)Nc1cc(C)nc2ccccc12</chem>
Expiration Date	12 months from date of receipt.