

Product Datasheet

kuwanon G (orb1303383)

Description	kuwanon G
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	C ₄₀ H ₃₆ O ₁₁
Purity	98.47%
MW	692.71
CAS Number	75629-19-5
SMILES	<chem>OC=1C(=C2C(C(=O)C(CC=C(C)C)=C(O2)C3=C(O)C=C(O)C=C3)=C(O)C1)[C@@H]4[C@@H](C(=O)C5=C(O)C=C(O)C=C5)[C@@H](CC(C)=C4)C6=C(O)C=C(O)C=C6</chem>
Expiration Date	12 months from date of receipt.