

Product Datasheet

Choline Fenofibrate (orb1303534)

Description	Choline Fenofibrate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{22}H_{28}ClNO_5$
Purity	99.95%
MW	421.91
CAS Number	856676-23-8
SMILES	<chem>CC(C)(C(=O)[O-])Oc1ccc(cc1)C(=O)c1ccc(cc1)Cl.C[N+](C)(C)CCO</chem>
Expiration Date	12 months from date of receipt.