

Product Datasheet

Protease-Activated Receptor-1, PAR-1 Agonist acetate (orb1303612)

Description	Protease-Activated Receptor-1, PAR-1 Agonist acetate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{37}H_{62}N_{10}O_{11}$
Purity	98.18%
MW	822.95
SMILES	<chem>CC(O)=O.O=C(N)C[C@@H](C(O)=O)NC([C@H](CCCNC(N)=N)NC([C@H](CC(C)C)NC([C@H](CC(C)C)NC([C@H](CC1=CC=CC=C1)NC([C@H]([C@H](O)C)N)=O)=O)=O)=O)=O</chem>
Expiration Date	12 months from date of receipt.