

Product Datasheet

Leukadherin-1 (orb1303971)

Description	Leukadherin-1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{22}H_{15}NO_4S_2$
Purity	98.49%
MW	421.49
CAS Number	344897-95-6
SMILES	<chem>OC(=O)C1=CC=C(C=C1)C1=CC=C(O1)\C=C1/SC(=S)N(CC2=CC=CC=C2)C1=O</chem>
Expiration Date	12 months from date of receipt.