

Product Datasheet

PF-8380 (orb1303978)

Description	PF-8380
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{22}H_{21}Cl_2N_3O_5$
Purity	99.45%
MW	478.33
CAS Number	1144035-53-9
SMILES	<chem>C(CCN1CCN(C(OCC2=CC(Cl)=CC(Cl)=C2)=O)CC1)(=O)C=3C=C4C(=CC3)NC(=O)O4</chem>
Expiration Date	12 months from date of receipt.