

Product Datasheet

Salubrinal (orb1304594)

Description	Salubrinal
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{17}Cl_3N_4OS$
Purity	99.97%
MW	479.81
CAS Number	405060-95-9
SMILES	<chem>C1C(Cl)(Cl)C(NC(=O)\C=C\C1=CC=CC=C1)NC(=S)NC1=CC=CC2=C1N=CC=C2</chem>
Expiration Date	12 months from date of receipt.