

Product Datasheet

Takeda103A (orb1305414)

Description	Takeda103A
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{23}F_2N_7O$
Purity	98.01%
MW	463.48
CAS Number	865609-72-9
SMILES	<chem>C(CC)N1C(=NN=C1CNC2=CC(C(NCC3=C(F)C=CC=C3F)=O)=CC=C2)C=4C=CN=CN4</chem>
Expiration Date	12 months from date of receipt.