

Product Datasheet

Tenofovir Disoproxil Fumarate (orb1305528)

Description	Tenofovir Disoproxil Fumarate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{19}H_{30}N_5O_{10}P \cdot C_4H_4O_4$
Purity	99.86%
MW	635.51
CAS Number	202138-50-9
SMILES	<chem>OC(=O)C=C\C(O)=O.CC(C)OC(=O)OCOP(=O)(CO[C@H](C)Cn1cnc2c(N)ncnc12)OCOC(=O)OC(C)C</chem>
Expiration Date	12 months from date of receipt.