

Product Datasheet

AdipoRon (orb1305879)

Description	AdipoRon
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{27}H_{28}N_2O_3$
Purity	99.63%
MW	428.52
CAS Number	924416-43-3
SMILES	<chem>C(=O)(C1=CC=C(OCC(NC2CCN(CC3=CC=CC=C3)CC2)=O)C=C1)C4=C C=CC=C4</chem>
Expiration Date	12 months from date of receipt.