

Product Datasheet

Tri-Salicylic Acid (orb1305916)

Description	Tri-Salicylic Acid
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{14}O_7$
Purity	99.14%
MW	378.33
CAS Number	85531-17-5
SMILES	<chem>O(C(=O)C1=C(O)C=CC=C1)C2=C(C(OC3=C(C(O)=O)C=CC=C3)=O)C=CC=C2</chem>
Expiration Date	12 months from date of receipt.