

Product Datasheet

Fipronil (orb1306142)

Description	Fipronil
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{12}H_4Cl_2F_6N_4OS$
Purity	98%
MW	437.15
CAS Number	120068-37-3
SMILES	<chem>NC=1N(N=C(C#N)C1S(C(F)(F)F)=O)C2=C(Cl)C=C(C(F)(F)F)C=C2Cl</chem>
Expiration Date	12 months from date of receipt.