

Product Datasheet

Telotristat (orb1306323)

Description	Telotristat
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{25}H_{22}ClF_3N_6O_3$
Purity	99.49%
MW	546.93
CAS Number	1033805-28-5
SMILES	<chem>[C@H](OC=1C=C(N=C(N)N1)C2=CC=C(C[C@@H](C(O)=O)N)C=C2)(C(F)(F)F)C3=C(C=C(Cl)C=C3)N4C=CC(C)=N4</chem>
Expiration Date	12 months from date of receipt.