

Product Datasheet

TGR5 Receptor Agonist (orb1307370)

Description	TGR5 Receptor Agonist
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{18}H_{14}Cl_2N_2O_2$
Purity	99.93%
MW	361.22
CAS Number	1197300-24-5
SMILES	<chem>C(N(C)C1=CC=C(Cl)C=C1)(=O)C=2C(=NOC2C)C3=C(Cl)C=CC=C3</chem>
Expiration Date	12 months from date of receipt.