

Product Datasheet

Tipifarnib (S enantiomer) (orb1307666)

| | |
|------------------------|--|
| Description | Tipifarnib (S enantiomer) |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{27}H_{22}Cl_2N_4O$ |
| Purity | 99.05% |
| MW | 489.4 |
| CAS Number | 192185-71-0 |
| SMILES | <chem>[C@@](N)(C=1C=C2C(=CC(=O)N(C)C2=CC1)C3=CC(Cl)=CC=C3)(C=4N(C)C=NC4)C5=CC=C(Cl)C=C5</chem> |
| Expiration Date | 12 months from date of receipt. |