

Product Datasheet

TAO Kinase inhibitor 1 (orb1307702)

Description	TAO Kinase inhibitor 1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{25}H_{24}N_2O_2$
Purity	99.45%
MW	384.47
CAS Number	850467-66-2
SMILES	<chem>N(C(CNC(=O)C1=CC=C(C=C1)C2=CC=CC=C2)=O)C3C=4C(CCC3)=CC=CC4</chem>
Expiration Date	12 months from date of receipt.