

Product Datasheet

SEC inhibitor KL-1 (orb1307748)

Description	SEC inhibitor KL-1
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{18}H_{16}ClNO_4$
Purity	100%
MW	345.78
CAS Number	900308-84-1
SMILES	<chem>N(C(C(=CC(=O)C1=CC(OC)=CC=C1)O)=O)C2=C(C)C=CC(Cl)=C2</chem>
Expiration Date	12 months from date of receipt.