

Product Datasheet

HI-TOPK-032 (orb1308229)

Description	HI-TOPK-032
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{20}H_{11}N_5OS$
Purity	98.69%
MW	369.4
CAS Number	487020-03-1
SMILES	<chem>C(#N)C=1C=2C(N3C1C=C(NC(=O)C4=CC=CS4)C=C3)=NC=5C(N2)=C=C5</chem>
Expiration Date	12 months from date of receipt.