

Product Datasheet

Dynarrestin (orb1308360)

Description	Dynarrestin
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{22}H_{23}F_2N_3O_2S$
Purity	99.49%
MW	431.5
CAS Number	2222768-84-3
SMILES	<chem>C(OCC)(=O)C1=C(N=C(NC2=C(F)C=C(F)C=C2)S1)C3=CC=C(N(CC)CC)C=C3</chem>
Expiration Date	12 months from date of receipt.