

Product Datasheet

Bemesetron (orb1308577)

Description	Bemesetron
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{15}H_{17}Cl_2NO_2$
Purity	97.37%
MW	314.21
CAS Number	40796-97-2
SMILES	<chem>O(C(=O)C1=CC(Cl)=CC(Cl)=C1)[C@H]2C[C@@]3(N(C)[C@](C2)(CC3)[H])[H]</chem>
Expiration Date	12 months from date of receipt.