

Product Datasheet

Atorvastatin lactone (orb1308630)

Description	Atorvastatin lactone
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{33}H_{33}FN_2O_4$
Purity	96.60%
MW	540.62
CAS Number	125995-03-1
SMILES	<chem>C(NC1=CC=CC=C1)(=O)C=2C(=C(N(CC[C@@H]3C[C@@H](O)CC(=O)O3)C2C(C)C)C4=CC=C(F)C=C4)C5=CC=CC=C5</chem>
Expiration Date	12 months from date of receipt.