

Product Datasheet

Tafenoquine Succinate (orb1309033)

Description	Tafenoquine Succinate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{28}H_{34}F_3N_3O_7$
Purity	99.84%
MW	581.58
CAS Number	106635-81-8
SMILES	<chem>O=C(O)CCC(O)=O.CC(NC1=C2N=C(OC)C=C(C)C2=C(OC3=CC=CC(C(F)(F)F)=C3)C(OC)=C1)CCCN</chem>
Expiration Date	12 months from date of receipt.