

Product Datasheet

PXS-5153A (orb1309203)

Description	PXS-5153A
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{20}H_{25}Cl_2FN_4O_2S$
Purity	99.71%
MW	475.41
CAS Number	2125956-82-1
SMILES	<chem>CC1=C(C2=CC(S(N(C)C)(=O)=O)=CC=C2)C3=NC=CC=C3N1C/C(F)=C/CN.Cl.Cl</chem>
Expiration Date	12 months from date of receipt.