

Product Datasheet

N6-(2-Phenylethyl)adenosine (orb1309362)

| | |
|------------------------|---|
| Description | N6-(2-Phenylethyl)adenosine |
| Storage | Store at -20°C for 12 months. |
| Note | For research use only |
| Formula | $C_{18}H_{21}N_5O_4$ |
| Purity | 99.94% |
| MW | 371.39 |
| CAS Number | 20125-39-7 |
| SMILES | <chem>N(CCC1=CC=CC=C1)C2=C3C(N(C=N3)[C@@H]4O[C@H](CO)[C@@H](O)[C@H]4O)=NC=N2</chem> |
| Expiration Date | 12 months from date of receipt. |