

Product Datasheet

ML303 (orb1309409)

Description	ML303
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{16}F_3N_3O_2$
Purity	99.79%
MW	399.37
CAS Number	1638211-04-7
SMILES	<chem>CC=1C=2C(N(N1)C3=CC=C(C(F)(F)F)C=C3)=NC=CC2C4=CC(OC)=C(O)C=C4</chem>
Expiration Date	12 months from date of receipt.