

Product Datasheet

Lisinopril dihydrate (orb1310366)

Description	Lisinopril dihydrate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{35}N_3O_7$
Purity	99.87%
MW	441.52
CAS Number	83915-83-7
SMILES	<chem>O.NCCCC[C@H](N[C@@H](CCc1cccc1)C(=O)O)C(=O)N1CCC[C@H]1C(=O)O.O</chem>
Expiration Date	12 months from date of receipt.