

Product Datasheet

Fenoprofen calcium dihydrate (orb1310858)

Description	Fenoprofen calcium dihydrate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{30}H_{30}CaO_8$
Purity	99.29%
MW	558.63
CAS Number	53746-45-5
SMILES	<chem>O.O.[Ca++].CC(C([O-])=O)c1cccc(Oc2ccccc2)c1.CC(C([O-])=O)c1cccc(Oc2ccccc2)c1</chem>
Expiration Date	12 months from date of receipt.