

Product Datasheet

Loteprednol etabonate (orb1310863)

Description	Loteprednol etabonate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{31}ClO_7$
Purity	99%
MW	466.95
CAS Number	82034-46-6
SMILES	<chem>O(C(OCC)=O)[C@@]1(C(OCCl)=O)[C@]2(C)[C@@](CC1)([C@]3([C@]([C@@H](O)C2)([C@]4(C)C(CC3)=CC(=O)C=C4)[H])[H])[H]</chem>
Expiration Date	12 months from date of receipt.