

Product Datasheet

Sorafenib (orb1310944)

Description	Sorafenib
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{21}H_{16}ClF_3N_4O_3$
Purity	99.61%
MW	464.82
CAS Number	284461-73-0
SMILES	<chem>O(C=1C=C(C(NC)=O)N=CC1)C2=CC=C(NC(NC3=CC(C(F)(F)F)=C(Cl)C=C3)=O)C=C2</chem>
Expiration Date	12 months from date of receipt.