

Product Datasheet

Sorafenib tosylate (orb1310945)

Description	Sorafenib tosylate
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{28}H_{24}ClF_3N_4O_6S$
Purity	99.81%
MW	637.03
CAS Number	475207-59-1
SMILES	<chem>Cc1ccc(cc1)S(O)(=O)=O.CNC(=O)c1cc(Oc2ccc(NC(=O)Nc3ccc(Cl)c(c3)C(F)(F)F)cc2)ccn1</chem>
Expiration Date	12 months from date of receipt.